

ENSURA™ High Purity Acids

Ensures to meet the demand of
ICP/ICPMS trace metal analysis
as per USP<232>, <233> & ICH Q3D

Product Feature :

- Available in 500 ml., 1 Ltr. & 2.5 Ltr. pack
- Suitable for all samples including Pharmaceutical, environment & food
- Best Suitable for sample preparation for ICP & ICPMS
- Stringent specification limit to make sure the consistent performance
- Tested for 68 elements in lower ppb level



Product Code	Product Description	Pack Size
20233-500ML	Hydrochloric acid ENSURA for trace metal analysis	500 ml.
20233-1LT	Hydrochloric acid ENSURA for trace metal analysis	1 Ltr.
20267-500ML	Nitric acid ENSURA for trace metal analysis	500 ml.
20267-1LT	Nitric acid ENSURA for trace metal analysis	1 Ltr.
20284-500ML	Sulphuric acid ENSURA for trace metal analysis	500 ml.
20284-1LT	Sulphuric acid ENSURA for trace metal analysis	1 Ltr.
20027-500ML	Hydrogen peroxide ENSURA for trace metal analysis	500 ml.